

HANDY

Recumbent/Upright

Owner manual



Handy Recumbent Assembly Procedures

Step 1:

Take out all parts in the box: (1) seat assembly (seat chasis, cushion pad, seat fabric support), (2) rear wheels x 2, (3) front frame assembly, (4) foot rest assembly, (5) main frame and (6) flag.

Accessories:

Upper left to lower right (7) seat belt, (8) fixing plate w/bolts for foot rest and rear axle, (9) foot rest velcro strap x 2, (10) shifter instruction, (11) rear wheel axle and hub cap x 2, (12) paint bottle, (13) steering spring set and (14) semicircle main frame fixing plate and clamp set.



Step 2:



Assemble the rear axle to the main frame.

Step 3:



Insert the skewer thru the rear hub and into the rear axle, put on the hub cover.

Step 4:



Assemble the hand pedals, gear shifter at right.

Step 5:



Assemble main frame to front frame assembly, clamp top and bottom.

Step 6:



Turn the front wheel left and insert the spring set, turn the wheel back and tighten.

Step 7:



Assemble the foot rest and apply the velcro strap.

Step 8:



Assemble the seat chasis to the main tube at the desired position and tighten all bolts.

Step 8 (cont'd):



Apply the seat fabric support, seat belt and cushion pads.

Step 9:



Set up the red flag and you are ready to go. Before riding, make sure all bolts are properly tightened.

Important Safety Checklist!

- Check and make sure that all clamps and QR are properly tightened.
- Always wear a helmet and abide by all traffic regulations and courtesies!

Important Maintenance Checklist!

- 1st service inspection after at most 300 km or at latest one months after purchase, repeat every 3 months.
- Parts need lubrication: chain, headset, bearings.
- Parts prohibit lubrication: seatpost, seat tube, adjustable stem, rims.
- Repeat lubrication every half month or after at most 150 km ride.
- Apply lubrication immediately after washing or riding in rainy day.